

Compatibility reaction tubes

The 48, 96, and 384 well blocks of the LabCycler are compatible to the most standard test tubes and microtiterplates. In 1996, the Society for Biomolecular Screening (SBS) began an initiative to create a standard definition of a microtiter plate. A series of standards was completed in 2003 and published by the American National Standards Institute (ANSI) on behalf of the SBS. The standard governs various characteristics of a microtiterplate including well dimensions (e.g. diameter, spacing and depth) as well as plate properties (e.g. dimensions and rigidity).

Compatibility 48 well block, 0.5 ml single tubes	
Abgene	AB-0350
Bio-Budget Technologies GmbH	my-Budget PCR-Tubes 0.5 ml, cat. no. 20-3-110-50-0
Greiner bio-one	Cat. No. 682201
Kisker Biotech	G-005 (serial)

Compatibility 96 well block, 0.2 ml single tubes	
Abgene	AB-0620
Bio-Budget Technologies GmbH	my-Budget PCR-Tubes 0.2 ml, cat. no. 20-3-010-50-0
Greiner bio-one	Cat. No. 683201
Kisker Biotech	G-001 (serial)

Compatibility 96 well block, 0.2 ml stripes	
Abgene	AB-0266, AB-0624, AB-0648, AB-0866, AB-1183
Bio-Budget Technologies GmbH	my-Budget 8er-Strips, cat. No. 30-SP-0020
Greiner bio-one	Cat. No. 673210
Kisker Biotech	G-002 (serial)

Compatibility 96 well block, microtiterplates	
Abgene	AB-0600, AB-0700, AB-0731, AB-0800, AB-900, AB-990, AB-1100, AB-1300, AB-1400, AB-1500, AB-1900 AB-2800
Bio-Budget Technologies GmbH	my-Budget "96 well" PCR-plates, cat. no. 30-SP-9600 my-Budget "96 well" PCR-plates, (Skirted), 30-SP-0037 my-Budget "96 well" PCR-plates, (Semi-Skirted), 30-SP-0127 my-Budget "96 well" PCR-plates, (Low Profile), 30-SP-0029
Greiner bio-one	Cat. No. 651501, 651550, 651560, 651570, 651585, 652201, 652270
Kisker Biotech	G-060 (serial)
Sorenson	Semi-Skirt 96 PCR plate Skirted 96 PCR plate 96-well ultra plate

Compatibility 384 well block, microtiterplates	
Abgene	AB-0937, AB-1111, AB-1384, TF-0384
Bio-Budget Technologies GmbH	my-Budget "384 well" PCR-plates, cat. no. 30-SP-0448
Greiner bio-one	Cat. No. 785201, 785290
Kisker Biotech	G-034
Sorenson	384-well PCR plate NX